

**PHCOG MAG:
Analysis of Plant Drugs
HPTLC Identification of Chaste
tree fruits (*Vitex agnus-castus*)
Courtesy : CAMAG**

Method: HPTLC

Sample for Identification: Hydrophilic flavonoids

Sample preparation:

0.75 g of powdered drug (fruits) is extracted for 2 min with 30 mL methanol in a turbo-extractor (polytron). After centrifuging for 5 min the supernatant is removed and set aside. The sediment is again extracted for 2 min with 30 mL methanol in the polytron, then centrifuged for 5 min. The combined supernatants are evaporated to dryness and reconstituted with 5 mL methanol. This is the test solution.

Standards (optional) preparation:

1 mg chlorogenic acid and 2 mg of iso-orientin are dissolved in 5 mL methanol.

Derivatization reagents:

Natural Products reagent (NP reagent): 1 g of diphenylborinic acid aminoethylester is dissolved in 200 mL ethyl acetate.

Macrogol reagent: 10 g of polyethylene glycol 400 (macrogol) are dissolved in 200 mL dichloromethane.

Chromatographic conditions:

Stationary phase:

HPTLC plates silica gel 60 F254 (Merck), 10x10 cm or 20x10 cm.

Mobile Phase:

Tetrahydrofurane, toluene, formic acid, water (8:4:1:0.5)

Sample application:

2 µL test solution and standard are applied as 8 mm bands, min. 2 mm apart, 8 mm from lower edge of plate.

Development:

CAMAG 10x10 cm or 20x10 cm Horizontal Developing Chamber (HDC), unsaturated, 2 mL (respectively 5 mL) developing solvent.

Developing distance 60 mm from lower edge of plate. The plate is then dried with a hair dryer (warm air) for 5 min.

Detection:


a) UV 254 nm

b) NP/Macrogol reagent: the plate is heated at 100°C for 3 min, then dipped while still hot in NP reagent, dried in a stream of cold air, then dipped in Macrogol reagent. Examination under UV 366 nm.

Results:

Chromatograms obtained are compared as shown below:

compare to the chromatograms below.



a) UV 254 nm

b) NP-reagent/Macrogol, UV 366 nm

Track 1: Chlorogenic acid,

Track 2: *Vitex agnus castus* (Croatia),

Track 3: iso-orientin

a) UV 254 nm

The reference solution shows a dark band for iso-orientin at $R_f = 0.31$ and a dark band for chlorogenic acid at $R_f = 0.41$. The test solution shows strong dark bands at $R_f = 0.31$ (isoorientin), $R_f = 0.36$, and at $R_f = 0.88$. There is a faint band at the position of chlorogenic acid and for additional weak bands between $R_f = 0.45$ and $R_f = 0.8$.

b) NP-reagent/Macrogol, UV 366 nm

The reference solution shows a strong yellow band for iso-orientin at $R_f = 0.31$ and a light blue band for chlorogenic acid at $R_f = 0.41$. The test solution shows bands matching those of iso-orientin and chlorogenic acid in color and R_f . There is a faint band at $R_f = 0.07$, a yellow brown band at $R_f = 0.45$, a olive green band at $R_f = 0.51$, a dark blue band at $R_f = 0.70$, a light blue band a $R_f = 0.78$ and an orange band at $R_f = 0.88$

NOTE: The flavonoid pattern of various samples of *Vitex agnus castus* does not vary significantly.

However, the content of flavonoids of some samples may be very low.

COURTESY : CAMAG APPLICATION NOTES

PHCOG MAG. is thankful to Dr. Marcel Hug, Camag, Switzerland for permission to use Application Notes.